



PACBIO®

# Experimental design and obtaining DNA for long-read sequencing



# Outline

- Select your species and/or individual
- Choose sequencing platform
- Generate high quality DNA
- Some examples





# What do all these things have in common?


Shifts in pollinators




Evolution of carnivorous plants



Response to nutrients





Flower color morphs






# Genome as a puzzle



# Genome size of angiosperms



*Arabidopsis thaliana*  
135 MB



*Oryza sativa*  
430 MB



*Zingiber officinale*  
1,582 MB



*Hordeum vulgare*  
5,300 MB



*Allium cepa*  
16,000 MB



*Tulipa sylvestris*  
59,241 MB



# Chromosome number and ploidy

**Kewscience** Royal Botanic Gardens Kew.org

## Plant DNA C-values Database

Home

### Plant DNA C-values Database

Release 7.1, April 2019. Leitch IJ, Johnston E, Pellicer J, Hidalgo O, Bennett MD  
<https://cvalues.science.kew.org/>

**All Plants** **Angiosperm** **Gymnosperm** **Pteridophyte** **Bryophyte** **Algae**

**VALUE**

The DNA amount in the unreplicated gametic nucleus of an organism is referred to as its C-value, irrespective of the ploidy level of the taxon. The Plant DNA C-values Database currently contains C-value data for 12,273 species comprising 10,770 angiosperms, 421 gymnosperms, 303 pteridophytes (246 ferns and fern allies and 57 lycophytes), 334 bryophytes, and 445 algae.

If you have comments and/or suggestions contact [dnac-value@kew.org](mailto:dnac-value@kew.org)

- Home
- Introduction
- Search
  - All Plant C-values
  - Angiosperm C-values
  - Gymnosperm C-values
  - Pteridophyte C-values
  - Bryophyte C-values
  - Algal C-values
- Release History
- Related Databases
- Contacts

**CCDB** CHROMOSOME COUNTS DATABASE

Enter a genus or genus and species

[Home](#) [About](#) [Browse](#) [Services](#) [Add new counts](#) [Contact](#)

Prof. Itay Mayrose Lab - Plant Evolution, bioinformatics, & comparative genomics

The **Chromosome Counts Database (CCDB, version 1.47)** is a comprehensive community resource for plant chromosome numbers. CCDB aims to combine existing data **resources** into an extensive central database that will be updated regularly by the community.

Users and researchers are encouraged to contribute to the accuracy and completeness of the data in CCDB by **submitting new counts**, or **reporting erroneous counts**.

To start browsing for chromosome numbers, use the **Browse** page, or the search box above.

Recommended citation:  
CCDB is built upon many individual data collection efforts. In addition to the main citation detailed below we recommend citing the major resources where the downloaded data were first assembled.  
Main citation: Rice et al. 2015. The Chromosome Counts Database (CCDB) – a community resource of plant chromosome numbers. *New Phytol.* 206(1): 19-26.

**Abstract.**

### Browse

- Flowering Plants  
Angiosperms
- Conifers, cycads and allies  
Gymnosperms
- Ferns and fern allies  
Pteridophytes  
(monilophytes and lycophytes)
- Mosses and liverworts  
Bryophytes


<http://ccdb.tau.ac.il>

# Genome size estimate

- For many nonmodel systems, there are no entries in the Kew database
  - Even if the genus is there, genome size can vary between species

## Kmer (estimates)


- Cleaned Illumina data
- Jellyfish paired with GenomeScope or RESPECT



<b>k = 19</b>	k-mer coverage	28.0
<b>property</b>	<b>min</b>	<b>max</b>
Heterozygosity (%)	3.64	3.65
Genome Haploid Length (bp)	11,995,570	12,010,675
Genome Repeat Length (bp)	2,179,917	2,182,662
Genome Unique Length (bp)	9,815,653	9,828,014
Model Fit (%)	98.26	98.89
Read Error Rate (%)	0.13	0.13


## Flow cytometry (more reliable)

- Fresh or silica dried material; protocols can vary a lot between species
- Need accurate references to compare



# Ideal scenario for sample selection



- Individual with lots of fresh material available
  - Single individual for sequencing and assembly
  - Scaffolding and/or annotation material can come from different individuals/species
  - Generate large amounts of high molecular weight DNA (often multiple micrograms)
  - RNA from multiple tissues and/or developmental stages
  - Fresh tissue for Hi-C sequencing
- "Clean" – less exposure to microorganisms or other organisms
  - If others are sequenced (including yourself), won't scaffold and be annotated





# Obtaining Sequencing data

Which puzzle is easier to put together?






# Which puzzle is easier to put together?




# In the world of genomes



Reference  
Genome


Sequencing  
Reads





# Short vs Long-reads

- Short reads
- Amplification errors and bias
- Several enzymatic steps
- Multi-molecule raw accuracy
- Errors tend to be systematic
- More coverage required



- Long reads
- No required amplification
- Simple sample prep
- Single molecule raw accuracy
- Errors tend to be random (vs. systematic)
- Less coverage required

# Long-read options


	PacBio Revio	SBS sequencing	Nanopore sequencing
Read length	15–20 kb	2x150 bp	10–100 kb
Read accuracy	99.95% (Q33)	99.92% (Q31)	99.26% (Q21)
Run time	24 hours <sup>3</sup>	44 hours	72 hours
Yield	90 Gb <sup>2,5</sup>	2,400–3,000 Gb	50–110 Gb
Variant calling – SNVs	✓	✓	✓
Variant calling – indels	✓	✓	✗
Variant calling – SVs	✓	✗	✓
5mC methylation	✓	✗	✓
Phasing	✓	✗	✓

HiFi targets 10 kbp, while Nanopore works “best” around 10-20 kbp (“best” can vary if fragment size or output is most desired)

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%



# Nanopore




# Nanopore approaches

## Genomic DNA

## RNA




DNA extraction  
(with incidental fragmentation)



Library Prep ~2.5 hours

Output 10-20 GB in 96 hours



10 minutes

8-10 GB in 96 hours



Full-length RNA

Primer annealing

Reverse transcription  
and strand switching

PCR with rapid  
attachment primers

Attachment of rapid 1D  
sequencing adapters

Loading

165 minutes

5 – 7 Million transcripts in 48 hours



# Getting started with Nanopore MinION



- Starter pack - \$1,999
  - Two 10.4.1 flow cells and 1 6-reaction library kit
- Flow cells – typically guaranteed to last 3 months or less in the fridge
  - Individually \$700 each
  - 24 bundle \$500 each
  - 48 bundle \$475
- Library prep is ~\$150 per sample (kit is \$100/sample + extra reagents)
  - 1.5-2 hour prep time, need 1 ug starting DNA
- PromethION is \$900 per flow cell, similar prep method, with increased sequencing output



# PacBio



SMRT Cell




Science, Vol 299, Jan 31 2003, pp682-686  
J. Appl. Phys. 103, 034301 (2008)

1. generate amplicon

2. ligate adaptors

3. sequence

4. data analysis



Fichot and Norman 2013; *Microbiome*



# Revio

- 90% of bases  $\geq$ Q30 and median read accuracy  $\geq$ Q30
- 15x increase in throughput over the Sequel II system
- Little less than 100 GB per SMRT cell
  - If DNA fragments are less than 10 KB, total output drops
- HiFi sequencing provides structural variants, repeat expansions, methylation, and haplotype phasing from a single library
  - **The \$1000 complete, phased genome**
- Typically outsourced as opposed to what can be done inhouse with minION




# Cost of sequencing genomes

- 50X Illumina:
    - 50Gb x \$26.5/Gb = **\$1,325**
  - 50X nanopore:
    - 50Gb x \$40/Gb = **\$2,000**
- 
- \$3,325**



# Hi-C

- Hi-C = high throughput chromatin conformation capture
- DNA of the same chromosome will be ***spatially*** close



# Scaffolding (1 GB genome)


- Hi-C Library kit \$500 + 50x Illumina \$550 = **\$1,050**
  - Library prep is not trivial, two-day protocol
  - Comes in sets of two
- Outsource to Phase Genomics
  - Send frozen samples on dry ice
  - Library prep \$1,500 + 150 M Read-pair Illumina \$750 = **\$2,250**
  - Guaranteed to get usable data
- Arima Genomics has a 6 hour rapid protocol
- Optical mapping by Bionano
  - Outsource (HWM extraction + Saphyr chip + analysis) = **~\$3,000**



# Extracting good DNA

# Tips for nonmodel systems

## Tissue Grinding




### Mortar and pestle

- Better yield, larger fragments
- Nanopore flow cell generated 18.5 GB with an N50 of 6.5 kb

### Grinding beads

Much lower yield; highly fragmented DNA  
Nanopore flow cell generated 12 GB with an N50 of 4.2 kb



## Extraction method

### SDS

(Monocots)



### CTAB



(Monocots)




Modified SDS for Monocots  
Modified CTAB for Eudicots

# Kit based approaches

- Several options available
- Zymo's is supposed to work in 45 minutes



DNA up to 150+ kb




## Quick-DNA HMW MagBead Kit

Cat #	Name	Size	Price <sup>1</sup>	Quantity
D6060	Quick-DNA HMW MagBead Kit	96 Preps	\$311.60	- 1 +




# QC

- Quantification – Need to rely on Qubit
  - Nanodrop drastically overestimates concentrations
  - 1 ug for each sequencing run; if size selection need around 5 ug starting out
- Purity/Cleanliness – Nanodrop
  - 260/280 values should be 1.8-2, while 260/230 values 2.0-2.2
  - If pure DNA, concentrations should be close to 1:1 (Nanodrop:Qubit)
- Integrity
  - Bioanalyzer or Femto Pulse
  - Low percentage agarose gel (0.5-1%) with low voltage
  - NEB 1 KB Extend Ladder (top band 48.5 kb)




# Size selection

- Blue pippin
- Circulomics Short Read Eliminator Kit
  - Target cutoff size: XS (10 kb), Regular (25 kb), XL (40 kb)
- Some/most DNA will be lost, but what remains is highly valuable
  - Be prepared for around a 40% quantity reduction with each cleaning step



# Cleaning the DNA



- Some species can be very difficult to get pure DNA
  - 1:1 Nanodrop:Qubit
- DNAeasy ProClean kit increases sequencing yield
  - DNA is sheared somewhat





# A couple of examples and costs

## Small(ish) genome

- Estimated genome size  
750 MB – 1 GB
- Two Nanopore flow  
cells (\$1,200)
- 50x Illumina (\$570)
- 1 SMRT cell Revio  
(\$2,760)
- 30x Hi-C (\$1,880)
- Total: **\$6,410**
- Chromosome scale with  
90% of estimated size in  
appropriate number of  
scaffolds

## Mid sized genome

- Estimated genome size  
1.8 GB
- Four Nanopore flow  
cells (\$2,400)
- 50x Illumina (\$570)
- 3 RSII SMRT cells  
(\$7,235)
- 30x Hi-C (\$2,250)
- Total: **\$12,545**
- Chromosome scale with  
90% of estimated size in  
appropriate number of  
scaffolds

## Larger genome

- Estimated genome size  
3 GB
- 4 RS II SMRT cells  
(\$9,090)
- 30x Hi-C (\$3,500)
- Total: **\$12,590**
- Chromosome scale with  
90% of estimated size in  
appropriate number of  
scaffolds

# Questions



@JLandisBotany



jbl256@cornell.edu